



## PCS 5000 Oxidizing Disinfectant Cleaning and Disinfection Process

### Features

#### PCS cotton cleaning cloths

- Cleans reduces surface soiling without the use of chemicals as effectively as cleaning with micro fiber cloths or cleaning with neutral or disinfectant cleaner.
- Use to clean, rinse, dry and apply disinfectants.
- Rapidly releases soils when rinsed or laundered.
- Cloths size and weight are large to increase soil pick up and are ideally suited for squirt bottle application of bleach disinfectants to dry cloths.
- 100 % organic

#### PCS 5000 Oxidizing Disinfectant/ Disinfectant Cleaner

- Equivalent to 10 to 1 bleach solution.
- Health Canada DIN: 02360500 registered
- Ready to use and stable formulation with a minimum of 24 months shelf life.
- Convenient flip top cap allows for easy application to PCS cotton cleaning cloths or surfaces like toilet bowls.

#### Benefits of the Process

- Reduces the formation of reactive by products of using bleach as PCS process removes organic soils and incompatible chemical residues prior to application of bleach solution.
- More effective cleaning than using pre moistened wipers.
- Eliminates mixing bleach.
- Safer than using 1 and 10 bleach solution in open bucket.
- More reliable than mixing 1 and 10 bleach solution and placing in work place bottles.
- Eliminates the issues of pre moistened wipers going dry prematurely as staff often times leave lids open.
- Re using cleaning and disinfectant application cloths is Environmentally Sustainable.
- Simple process with minimal staff training required
- Cost effective and substantially less expensive than using pre moistened wipers.
- Safer process for improved staff acceptance to the use of 10 to 1 bleach solutions,
- Easily demonstrated and validated with ATP testing.



# YOU WILL NEED

## PCS 5000 Oxidizing Disinfectant 946 ml container Laundered PCS Cotton Cleaning Cloths Acrylic bowl mop and caddy Cleaning bucket with warm water

### The Process:

- 1) Remove visible soils by applying PCS Friction and wiping dry with PCS cotton cloths  
It is important to remove urine, soap films and other organic soils before applying any disinfecting solution.
- 2) Squirt PCS 5000 Disinfectant/Disinfectant Cleaner onto PCS cotton cloth then wipe surfaces to clean/disinfect. Allow dwell time. Rinse or allow to air dry.
- 3) Clean toilet bowls with PCS 5000 Disinfectant/Disinfectant Cleaner.  
Use with Acrylic bowl mops.

Follow your institutions policy and procedures.

Insure laundry processes are effective in removing soil and disinfecting cotton cloths, process doe's not include addition of fabric softeners, anti chlor agents and that sufficient rinse cycles are used to insure cloths are clean and free of chemical residues.

Read material safety data sheet before use .Always wear recommended personal protective equipment. Not all surfaces are compatible with 5000 parts per million of bleach.

**5955-6 PCS 5000 Oxidizing  
Disinfectant/Disinfectant Cleaner  
6 x 946 ml**

**PCS Cotton Cleaning Cloths  
40 cm x 48 cm 60 per package**



Affix Business Card Here

**PCS**  
PROCESS CLEANING SOLUTIONS  
Bringing Science to the Art of Cleaning™

[www.processcleaningsolutions.com](http://www.processcleaningsolutions.com)  
Toll Free: 877.745.7277