

PCS Green Cleaning and Disinfecting Processes to Reduce Hospital Acquired Infections

PCS happy hands program for improved hand hygiene compliance. ECOgent Moisturizing Hand Wash with Vitamin E.

- Soap Free Hand Washing Ideally suited for frequent hand washing. Soap and surfactant based ingredients in hand soaps can with frequent use contribute to skin dryness and irritation.
- ECOgent is applied to dry hands and rapidly and thoroughly covers all parts of hands and rinses off hands more rapidly than soap.
- Hand soaps contain soap and surfactants which contribute to the pollution of our waterways. ECOgent contains no soaps or surfactants.
- Independently tested and validated to perform on par with national brand hand soap.

5981-12 500 mL w pump x 12/cs 5981-4 3.78L x 4/cs 5981-6 2L w pump x 6/cs

ECOgent Waterless Hand Washing - Hand Washing Without Water.

- Apply to hands (Minimum of 3 mL)
- Rub in as you would soap (10 to 15 seconds)
- Wipe completely dry with soft paper towel (Take extra time to insure to absorb all liquid from hands.)
- Designed for use where access to potable water is limited.

5982-12 500 mL w pump x 12/cs 5982-6 2L w pump x 6/cs

PCS cleaning process for a cleaner health care environment.

Apply PCS Friction to surfaces and wipe dry with a PCS cotton cloth.

PCS Cotton Cleaning Cloths

- Use to clean, rinse, dry and apply disinfectants, more effective than microfiber cloths.

PCS Friction

- Natural Organic Multi Purpose Cleaner.
- Non hazardous cleaning certified for reduced environmental impact and suitability for use around chemically sensitive individuals.

Easy to use.

- Apply and wipe dry with PCS cotton cloths.
- Rapidly dissolves and removes daily water deposits, soap scum, organic soils including fecal matter, urine and most soils spread from daily human activity.
- Independently tested and validated to clean to the scientific standard proposed by the World Federation of Building Service Contractors.

5978-4 Concentrate closed loop dispenser 3.78L x 4/cs 5978-6 Concentrate bottle 946ml x 6/cs



"Cleaning is always essential prior to disinfection or sterilization. An item that has not been cleaned cannot be assuredly disinfected or sterilized."
Public Health Agency of Canada/Health Canada

PCS MicroClean

- MicroClean dissolves soils and enhances the friction between cleaning cloth and surfaces to remove an amazing amount of soil.
- Cleans a wide variety of soils found on floors, carpets, washrooms, glass tables walls furnishings even delicate equipment safely and effectively.
- Economical with a single dilution for all cleaning applications 1 part cleaner to 256 parts water.
- Natural cleaner containing buffered lactic acid, sodium citrate salt of citric acid and table salt.
- Safe and non hazardous by WHIMIS criteria.
- Certified for use around sensitive individuals and for reduced environmental impacts.

5943-6 Concentrated 946 ml bottle x 6/cs
5943-4 Open Stock 3.78 L x 4/cs

5942-4 Concentrate closed loop dispenser 3.78L x 4/cs
5940-6 Concentrate closed loop dispenser 2L x 6/cs



Cleaning to protect public health.

- Reduces the potential of spreading germs from cleaning as compared to cleaning with a damp cloth with diluted neutral, general purpose or disinfectant cleaner and leaving surfaces damp.

Apply Disinfectant to PCS cotton cloth then apply to pre cleaned surfaces where required.

PCS 250 Oxidizing Disinfectant/Disinfectant Cleaner. Guarantee .

0.025% sodium hypochlorite /250 ppm of available chlorine.

- PCS preferred disinfectant.
- Ready to use product in sealed 946 ml bottle with flip top.
- Stable with a shelf life of twenty four months.
- Eliminates mixing of disinfectants with dispensers.
- Safer process improving staff acceptance.
- PCS cleaning and disinfecting process provides more effective cleaning and
- Disinfecting of the health care environment.

5908-6 6x 946ml Health Canada DIN 02314843.



PCS 1000 Oxidizing Disinfectant/Disinfectant Cleaner. Guarantee .

0.1% sodium hypochlorite /1000 ppm of available chlorine.

- PCS preferred disinfectant.
- Ready to use product in sealed 946 ml bottle with flip top.
- Stable with a shelf life of twenty four months.
- Eliminates mixing of disinfectants with dispensers.
- Safer process improving staff acceptance.
- PCS cleaning and disinfecting process provides more effective cleaning and
- Disinfecting of the health care environment.

5906-6 6x 946ml Health Canada DIN 02314851.



PCS 5000 Oxidizing Disinfectant / Disinfectant Cleaner.

Guarantee .5 % sodium hypochlorite /5000 ppm of available chlorine.

- Ready to use product in sealed 946 ml bottle with flip top.
- Stable with a shelf life of twenty four months.
- Eliminates mixing of disinfectants with dispensers.
- Safer process improving staff acceptance.
- PCS cleaning and disinfecting process provides more effective cleaning and disinfecting of the health care environment.

5955-6 6 x 946 ml containers Health Canada DIN 0236050.



www.processcleaningsolutions.com
Toll Free: 877.745.7277