


**Transform Your Clean:
 PCS Buffered pH Bleach
 The Ultimate Oxidizing
 Disinfectant Cleaner!**



Low Odour, 1-Minute Oxidizing Cleaner
 Hospital-Grade Disinfectant
 Broad Spectrum Virucide & Fungicide
 Eliminates Clostridioides Difficile Spores in 5 Minutes

PCS Buffered pH Bleach Oxidizing Disinfectant Cleaner • DIN: 02548534
 Active ingredient Sodium hypochlorite 0.26 % w/w when packed

PCS Buffered pH Bleach: The Safer, More Effective Disinfectant

Surface and Staff Safe: With a pH of 8.5 and 50% less sodium hypochlorite, it maintains disinfecting power while being safer for surfaces, staff, and occupants. Kills Clostridioides difficile spores.

Versatile Application: Suitable for various application techniques, making it perfect for protecting public health through oxidizing cleaning and disinfecting.

Wide Range of Uses: Ideal for hard, non-porous surfaces in homes, food processing areas, hospitals, healthcare facilities, institutions, schools, animal housing, and hospitality industries to end and prevent pathogen outbreaks.

Critical Care Decontamination: Perfect for professional decontamination of food contact surfaces, washrooms, and animal husbandry, including veterinary clinics, and public settings at risk of pathogen transmission.

Biofilm Penetration: More effective than alkali bleach at penetrating and dissolving biofilm matrices, making it ideal for healthcare environments.

Prevent and Eliminate Biofilms: Regular use prevents and removes biofilms. Ideal for deep cleaning to end pathogen outbreaks and reduce pathogen transfer in healthcare settings.

Proudly Canadian: Designed and manufactured with advanced Canadian technology.

Choose PCS Buffered pH Bleach Oxidizing Disinfectant Cleaner for a safer, more efficient way to maintain a clean and healthy environment.

Bacteria, Virus, Fungus and C. Diff Spores	Contact time
Staphylococcus aureus (ATCC 6538)	1 minute
Pseudomonas aeruginosa (ATCC 15442)	1 minute
Feline Calicivirus (as surrogate for Norovirus)	1 minute
Adenovirus Type 5 (ATCC #VR 1516)	1 minute
Trichophyton interdigitale (ATCC #9533)	1 minute
Clostridioides difficile Spores (ATCC #43598)	5 minutes

Product Code

#6283-6 6 x 946 mL per carton
 #6283-4 4 x 3.78 L carton
 #6285-2 PCS Buffered pH Bleach Wipe Application Kit
 (contains 2 x 946 ml PCS Buffered pH Bleach
 2 only rolls of 80 wipes 25x 30.48 cm and 2 dispensing buckets)

Keep out of reach of children.
 Caution may irritate eyes.