



BUY CANADIAN



PCS 1000 Plus Oxidizing Disinfectant Cleaner Concentrate. Diluted solution are effective Hospital Grade Disinfectant and Broad spectrum virucide in 1 minute.

Save up to 90% versus buying Ready to Use disinfectants.

PCS 1000 Plus Oxidizing Disinfectant Cleaner Concentrate

Active Ingredient 2% w/w Sodium Hypochlorite when packed.

Oxidizing cleaner
Oxidizing hospital grade Disinfectant
Oxidizing broad spectrum virucide.

DIN: 02521504

When diluted, this product contains 0.13% sodium hypochlorite and 0.01% hypochlorous acid. Regular use of high-level chlorine test strips (to show 1400 ppm available chlorine) is recommended.

To disinfect frequently touched surfaces such as non-critical medical equipment, bed rails, washroom fixtures, and surfaces during outbreaks in healthcare facilities, long-term care settings, schools, and other institutions: Remove visible soil, then apply the diluted solution to the surfaces. Allow the surface to remain wet for **1 minute**. Wipe the surface dry, rinse, or allow to air dry.

