

PCS Buffered pH Bleach Oxidizing Disinfectant Cleaner



Low odour 1 minute oxidizing cleaner
 Oxidizing hospital grade disinfectant
 Oxidizing broad spectrum virucide and fungicide
 Kills spores of *Clostridioides difficile* in 5 minutes

DIN: 02548534

Active ingredient Sodium hypochlorite 0.26 % w/w when packed

- Safer for surfaces, staff, and building occupants, contains a pH of 8.5 and 50% less sodium hypochlorite while maintaining disinfecting claims; kills *Clostridioides difficile* spores.
- Can be safely used with multiple application techniques, making it an ideal choice for oxidizing cleaning and disinfecting to protect public health.
- PCS Buffered pH Bleach Oxidizing Disinfectant Cleaner for use on hard non-porous environmental surfaces in domestic settings, food processing areas, hospitals, healthcare facilities, institutions, schools, animal housing areas, and hospitality industries.
- Ideal choice for critical care decontamination, professional service decontamination of food contact surfaces and washrooms, decontamination in animal husbandry including veterinary clinics, and most public settings that may pose a risk of transmitting pathogens from fomites.
- Ideal choice for decontamination of healthcare environments. PCS Buffered pH Bleach Oxidizing Disinfectant Cleaner penetrates and dissolves biofilm matrix more effectively than alkali bleach.
- Use regularly to prevent and remove biofilms. Use for deep cleaning to prevent and end pathogen outbreaks. Use regularly where biofilms form in patient sinks; spray in drains. Use for discharge cleans to reduce pathogen transfer to the next patient occupying the room.
- Canadian technology designed and manufactured in Canada.

Bacteria, Virus, Fungus and C. Diff Spores	Contact time
Staphylococcus aureus (ATCC 6538)	1 minute
Pseudomonas aeruginosa (ATCC 15442)	1 minute
Feline Calicivirus (as surrogate for Norovirus)	1 minute
Adenovirus Type 5 (ATCC #VR 1516)	1 minute
Trichophyton interdigitale (ATCC #9533)	1 minute
Clostridioides difficile Spores (ATCC #43598)	5 minutes

Product Code

#6283-6 6 x 946 mL per carton
 #6283-4 4 x 3.78 L carton
 #6285-2 PCS Buffered pH Bleach Wipe Application Kit
 (contains 2 x 946 ml PCS Buffered pH Bleach
 2 only rolls of 80 wipes 25x 30.48 cm and 2 dispensing buckets)

Keep out of reach of children.
 Caution may irritate eyes.